Push® is a registered trademark of Net International bv

Medical Device

04/2020

47207 / 2

Made in the EU

Designed in the Netherlands

www.push.eu

The Netherlands

Eurospineen 31
Eindhoven 5600 AB
Netherlands

4.70.1

REF

PUSH SPORTS

ELBOW BRACE

Useful information for you!

PUSH SPORTS ELBOW BRACE

The Push Sports Elbow Brace was developed to address strains at the insertion points of the extensor and flexor muscle tendons in the lower arm. These strains (commonly referred to as tennis or golfer's elbow) are caused by a (sudden) sharp strain, or prolonged strain, on the muscles involved. The elbow brace is easy to apply, adaptable to your individual needs and has specialized features to ensure it won't slide during movement and activity.

WHEN TO USE THE PUSH SPORTS ELBOW BRACE?

The Push Sports Elbow Brace is used primarily to relieve the strain on tendons during high-strain activities, such as exercising and household or professional activities. The brace is proven to relieve tennis or golfer's elbow.

INDICATIONS

- Tennis elbow (lateral epicondylitis)
- Golfer's elbow (medial epicondylitis)

CONTRAINDICATIONS

- There are no known disorders where the Push Sports Elbow Brace should not be worn.

CAUTIONS

- Ensure that the brace is not fitted too tightly. Use the brace only for the appropriate indications. Proper use of your Push Sports brace may require diagnosis by a doctor or therapist. Pain is always a warning sign. If pain or complaints persist without a clear cause, it is recommended that you visit your GP, physiotherapist or sports physician.
- Push Sports braces were designed to be worn on an undamaged skin. Damaged skin must first be covered with suitable material.
- Optimal functioning can only be achieved by selecting the right size. Therefore, preferably fit the brace before purchasing, to check its size. In the case of an online purchase, always ensure that there is a good return service. Your reseller's professional staff will explain the best way to fit your Push Sports brace. We recommend that you keep this product information in a safe place, to be able to consult it for fitting instructions whenever needed.
- Check your Push Sports brace each time before use for signs of wear or ageing of components or seams. Optimal performance of your Push Sports brace can only be guaranteed if the product is fully intact.
- Within the framework of EU regulation 2017/745 on medical devices, the user must report any serious incident involving the product to the manufacturer and to the competent authority of the Member State in which the user is established.

SPECIAL FEATURES

The elbow brace uses a pad to exert pressure on the affected lower arm muscles and thus relieves the strain on the tendon insertions of these muscles to the elbow. The effectiveness of this brace is a function of the viscoelastic material used in the pad, and a small elastic section on the strap. These features cause the pressure to increase gradually during flexion or extension of the muscles.

TAKING MEASUREMENTS

The elbow brace is available in one size, universally applicable for the left and right sides.

FITTING

Before use: remove the protective foil from the pad.

- Insert the end of the strap through the strap ring. Slide the brace up the arm. In the case of a tennis elbow, the yellow pad is positioned on the outside of the elbow; in the case of a golfer's elbow, it should be on the inside. The strap ring should not be located in the crease of the elbow.
- Place the brace on the lower arm muscles, approximately six centimetres below the elbow.
- Close the strap tightly enough to feel pressure when the arm is bent and the muscles are tightened. No tension should be felt when the muscles are relaxed. The strap can be left in the strap ring when the brace is removed.

pushsports.eu/elbow-brace/applying

For details of the distributor in your country, we refer to www.push.eu/distributors

WASHING INSTRUCTIONS/MAINTENANCE

To promote hygiene and to maximise the lifespan of the brace, it is important to wash the brace regularly, observing the washing instructions. The brace can be washed at temperatures of up to thirty degrees Celsius, on a delicate fabric cycle or by hand. Before washing, fasten the hook and loop material to protect the brace and other laundry. Do not use any bleaching agents. Let the brace dry in the open air (**not** near a source of heat or in a tumble drier). Never store the brace when it is wet. Always let it dry first.

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