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Instructions for use

Carpalorthese

Carpal Orthosis

dynamics[®]

ofa bamberg

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EN english	Carpal Orthosis	Fitting by experts	6
		Important information	10
		Garantie	20

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Ofa Bamberg. **More than 90 years' experience** **with your health.**

Ofa Bamberg is one of the leading manufacturers of medical supplies in Germany. Since the company was founded in 1928, we have set extremely high standards for our products: maximum comfort and optimal efficacy.

In addition to medical supports and orthoses, Ofa Bamberg also offers an extensive range of medical compression stockings and preventive travel and support stockings.

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Moulding the plastic splint (by experts only)

The plastic splint is thermoformable and can be moulded individually for you by trained experts. This requires the inner padding and Velcro straps to be removed. The areas to be moulded are heated evenly to approximately 150 °C using a hot air blower. The splint can then be moulded into the desired shape then left to cool. The material can be heated back up and moulded again if necessary.

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Introduction

Dear user, we are delighted that you have chosen an Ofa Bamberg product and we want to thank you for putting your trust in us. Please take a moment to read this product information through carefully to ensure your Dynamics product works as effectively as possible.

Ofa Bamberg wishes you a speedy recovery!

Fitting instructions

The pictures that go with the fitting instructions can be found on the inside cover (2) at the start of the brochure.

1

To put the orthosis on, open all the Velcro fasteners and place your lower arm into the splint. Move your fingers and thumb carefully until your hand is in an optimal position.

2

Then close the Velcro strap around the wrist followed by the Velcro straps around the fingers and thumb and finally the Velcro strap for the lower arm. Ensure that the product is secure but not too tight so as to prevent pinching.

3

You can also put the additional strap for the back of the hand on. The padding should cover the sensitive back of the hand. Wrap the Velcro strap around the orthosis and close it by inserting it through the eyelet provided and fastening it tightly.

Intended Use

The Dynamics Carpal Orthosis

stabilises and immobilises your entire hand. The thermoplastic material restricts the movement of your hand and immobilises it in a functional position. This supports the healing process.

Important information

- The initial fitting and adjustment of the product as well as instructions on its proper and safe use should be provided by trained medical/orthopaedic specialists. → See page 6.
- The product should only be used for the indications in the adjacent column.
- Ensure that the product fits correctly and check it before each use for functionality, wear and damage (e.g. breaks, tears or deformations). A damaged product may no longer be worn.
- Wear the product directly on the skin. It should only be worn over intact skin. If it is to be worn in combination with other products, please discuss this with your doctor first.
- Unless otherwise prescribed, the product should be removed during long periods of rest (e.g. when sleeping).
- If you experience an unpleasant sensation, severe pain or other symptoms while wearing the product, remove it and contact your doctor or specialist supplier.
- We use a comprehensive quality assurance system to check our products. If you have any complaints, however, please contact your specialist supplier.

Instructions for care

- The product should be washed before its first use and then on a regular basis.
- Remove the detachable padding and Velcro from the splint to wash. Please ensure that the adhesive strips on the splint are not also removed. The splint can be cleaned with a damp cloth.
- Wash the textile parts by hand at a maximum of 30 °C. Use mild soapy water and rinse the product well. Squeeze the water out carefully without wringing the product.
- Mould the product into the right shape and let it dry in the air.
- To maintain the quality of your product, please do not use any chemical cleaning products, bleaches, petrol or fabric softener on it. Do not apply any fatty or acidic creams, ointments or lotions to the skin in the area where the product is used. These substances may damage the material.
- Store the product in a dry place and protect it from sunlight and heat, preferably in the original packaging.

Reporting obligation

In accordance with legal provisions within the EU, patients and users are obliged to report any serious incident that occurs when using a medical device to both the manufacturer and the competent national authorities (in Germany BfArM, Bundesinstitut für Arzneimittel und Medizinprodukte, the Federal Institute for Drugs and Medical Devices) immediately.

Indications

- Severe sprains and bruises
- Severe tendon sheath inflammation (tendinitis)
- Carpal tunnel syndrome (nerve damage in the hand area)
- Aftercare of breaks, joint injuries (dislocations) and operations
- Paralysis after a stroke (hemiplegia)


Side effects

There are currently no known side effects for any part of the body if the product is used correctly. However, if the product is too tight it can cause local pressure points or constrict blood vessels and nerves.

Contraindications

- Consult your doctor before using this product in the following cases:
- Skin diseases or injuries in the area where the product is used, particularly if you experience any signs of inflammation such as redness, an increase in temperature or swelling
 - Feeling of numbness and circulatory problems in the area where the product is used
 - Problems with lymphatic drainage and unexplained swelling

Care

-  Hand wash
-  Do not bleach
-  Do not tumble dry
-  Do not iron
-  Do not dry clean

Do not use fabric softener!

Material composition

40 % Velour
39 % PU
19 % PVC
2 % Nylon
POM

Disposal

Please dispose of the product according to local regulations after the end of use.

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Guarantee

If used as intended and maintained property, the product can be used for six months. The product should be handled carefully in line with these instructions for use. Improper use or improper changes to the product may negatively impact the performance, safety and functionality of the product and lead to the exclusion of claims. The product is designed for the treatment of a single patient. It is not designed to be reused by a second patient or any further patients.

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Ofa Bamberg GmbH