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**Flachstrick**  
Flat-knitted

the medical  
compression  
specialist  
**since**  
**1928**  
made  
in Germany

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The flat knitted  
medical compression  
with merino wool

**ofa<sup>®</sup>**

---

**Ofa.**

Competence in medical compression.  
Traditionally innovative.

**Ofa is one of the leading German manufacturers of medical compression stockings as well as supports and braces.**

And not without reason. For since the company was founded in the year 1928, Ofa has been totally clear about the demands it places on its products: maximum comfort with optimal effectiveness.

Apart from a large choice of compression stockings, Ofa also offers a comprehensive range of supports and braces (Dynamics and Dynamics Plus). Preventive health products, such as the support and travel stockings of the Gilofa brand, round off the product portfolio.

|                |                       |    |                                       |    |
|----------------|-----------------------|----|---------------------------------------|----|
| DE deutsch     | Allgemeine Hinweise   | 6  | Anziehtipps                           | 7  |
| EN english     | General information   | 8  | Tips for putting the stockings on     | 9  |
| NL niederlands | Algemeen advies       | 10 | Aantrektips                           | 11 |
| IT italiano    | Indicazioni generali  | 12 | Consigli per indossare le calze       | 13 |
| PT português   | Instruções gerais     | 14 | Dicas para calçar as meias            | 15 |
| EL ελληνικά    | Γενικές υποδείξεις    | 16 | Συμβουλές πώς να φορέσετε τις κάλτσες | 17 |
| TR türkçe      | Genel bilgiler        | 18 | Giyme önerileri                       | 19 |
| DA dansk       | Generelle anvisninger | 20 | Sådan tages strømperne på             | 21 |
| FI suomi       | Yleisiä ohjeita       | 22 | Pukemisivinkkejä                      | 23 |
| RU русский     | Общая информация      | 24 | Советы по надеванию                   | 25 |

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## **Lastofa Forte**

### **Our all-rounder in flat-knit**

Lastofa Forte, the flat-knitted medical compression with merino wool. It attracts with its unique material composition. In addition to the synthetic fibres essential for medically correct compression, our Multitalent consists of a large extent of the finest merino wool. Here is an overview of the most important advantages of Lastofa Forte:

- exceptional comfort
- regulate heat and moisture
- breathable, with odour control
- strong but yet very soft
- easier to put on and take off

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## Dear User,

we are delighted that you have chosen Ofa compression stockings and thank you for putting your trust in us. Please take a moment to read through this product information carefully. For optimal compression therapy effectiveness!

## Intended use

Compression stockings are used to treat phlebological and lymphological arm and leg conditions, as well as lipoedema. Different products are used depending on the indication.

## Effective principle of medical compression stockings

Products tailored to the respective indication and controlled pressure accelerate blood reflux, improve microcirculation in the tissue and increase lymph flow.

## Important information

Only wear compression stockings after consultation with your doctor. Wearing your stockings daily and holding regular consultations with your doctor will ensure successful treatment! The stockings must be obtained from a medically trained specialist retailer. Your specialist retailer will also advise you on how to put the compression stockings on correctly. If you experience leg pain or skin irritation while wearing these stockings, please consult your doctor immediately. You can prevent skin reactions in the area around the grip tops by slightly adjusting the grip tops several times during the day. Please contact your doctor about wearing the stockings in combination with using other medical devices. We test our products within the scope of a comprehensive quality assurance. However, if you still have cause for complaint, please contact your specialist retailer. Do not attempt to perform repairs yourself. Return stockings to the specialist retailer who sold them to you. Any pilling of the material that may arise due to friction with other garments will not reduce the stocking's quality.

## Care instructions

The product is suitable for immediate use. We recommend washing the product before wearing it for the first time. Then wash daily by hand or on a delicate fabric program. Wash the stockings separately with mild detergent without fabric softener (e.g. Ofa Clean special detergent) at a maximum temperature of 40 °C and rinse them well. The stockings can be spun gently, but do not wring them out. Dry your compression stockings in air or in a dryer on a gentle cycle. Do not place them on a heater or in direct sunlight. Never iron the stockings. To maintain the effectiveness of your stockings, please do not use chemical cleaning agents, petrol, or fabric softeners. These substances damage the knitted fabric. Compression stockings are largely resistant to fats and oils, but should not come into contact with ointments. Stockings with silicone-coated grip tops should generally not come into contact with solvents, ointments, re-greasing washing lotions and body creams. Otherwise, the adhesive strength of the grip tops could be

### Indications (see also AWMF Medical Compression Therapy Guideline 037/005, Version 12/18)

- Improvement of venous symptoms and quality of life with chronic venous diseases
- Prevention and treatment of venous oedema
- Prevention and treatment of venous skin alterations
- Eczema and pigmentation
- Lipodermatosclerosis and atrophie blanche
- Therapy and pain reduction for venous leg ulcers, including relapse prevention
- Treatment of arterial and venous leg ulcers (see contraindications!)
- Varicose veins
- Initial phase after treatment for varicose veins
- Functional venous insufficiency
- Venous malformations
- Venous thrombosis
- Condition after thrombosis
- Post-thrombotic syndrome
- Thrombosis prophylaxis in mobile patients
- Lymphoedema
- Oedema in pregnancy
- Post-traumatic/post-operative oedema
- Cyclic idiopathic oedema
- Stage II lipoedema and above
- Venous stasis resulting from immobility (arthrogenic congestive syndrome, paresis, and partial paresis of the extremity)
- Occupation-related oedema (occupations involving long periods of standing and sitting)

- Medication-related oedema, if change of medication not possible
- Obesity with functional venous insufficiency
- Inflammatory dermatoses of the legs
- Nausea, dizziness, and congestion complaints during pregnancy
- Condition after burns
- Scar treatment

### Contraindications

#### Absolute contraindications:

- Advanced peripheral arterial occlusive disease (if any of the following values apply: ABPI < 0.5, ankle arterial pressure < 60 mmHg, toe pressure < 30 mmHg or TcPO<sub>2</sub> < 20 mmHg on instep)
- Congestive heart failure
- Septic phlebitis
- Phlegmasia cerulea dolens

#### Risks to be considered (therapy decision to be made based on benefit/risk assessment)

- Weeping dermatoses
- Allergies to compression stocking materials
- Severe sensitivity disorders of the extremity
- Advanced peripheral neuropathy (e.g. diabetes mellitus)
- Chronic polyarthritis

**No liability is accepted if these contraindications are not observed.**

## Risks and side effects

Medical compression stockings can cause skin necrosis and pressure damage to the peripheral nerves. Padding is recommended in at-risk areas. Users with sensitive skin may experience itching and skin irritation when using compression stockings. Adequate skin care is therefore recommended.

If you suffer the following symptoms, please remove the stockings immediately and discuss further treatment with your doctor: blue or white coloration of the toes, sensations of discomfort and numbness, increasing pain, shortness of breath and sweating, acute movement restrictions.

### Warranty

The shelf life of medical compression stockings is limited to 3 years when properly stored. Added to this is the wearing time, which is a maximum of 6 months with regular use and proper care. After this period, the products are excluded from the warranty. Defects attributable to normal wear and tear or improper use are not covered by the guarantee. Improper use also includes any subsequent modification to the product that has not been carried out by the manufacturer.

Slight colour deviations, which may arise due to the individual dyeing process, are not considered to be a defect. The product is intended for use by a single patient. If it is passed on to other persons for reuse, the manufacturer's product liability is no longer valid.

### Disposal

After the end of the usage period, please dispose of the product in accordance with local regulations.

### Mandatory reporting

Due to legal requirements, patients and users are obliged to report immediately any serious incidents which have occurred in connection with the product and which have resulted in, or could have led to a significant deterioration in the state of health or to death, to the manufacturer and the competent national authority (in Germany BfArM).

## english

permanently impaired (leg hair and dry skin can also limit the function of the grip tops). Store the stockings in a dry place away from sun and heat, in the original packaging if possible.

**Please wash right way out and do not use fabric softener!**

-  Wash on a delicate or gentle cycle
-  Do not bleach
-  Dry gently and carefully
-  Do not iron
-  No chemical dry cleaning


**Material composition**

Polyamide (Polyamid), elastane (Elastan) and wool (Wolle) are used. Please refer to the textile label sewn into your stocking for the precise details.


**Tips for putting the stockings on**

Put on the stockings first thing in the morning after getting up, as long as your legs are free from swelling. Please do not pull the compression stockings up by the edge or the grip top, instead, slide it upwards with your hands. We recommend wearing household rubber gloves to make it easier to put on the stockings and to protect the knitted fabric. Jewellery and long fingernails can damage the stockings. Please also pay attention to regular foot care and proper footwear. Ask your medical retailer about Ofa Grip special gloves and the practical donning and doffing aids for open and closed toe compression stockings.


**How to put them on, closed toe**

To put the stockings on, turn it inside out as far as the foot. Slip the foot section over your toes and then over the instep.

Now carefully pull the leg section over your heel. Please avoid pulling up the stockings at the end or grip top.


Pull the rest of the stockings up evenly without overstretching it. Finally, adjust the stockings on the leg until there are no wrinkles.

**How to put them on, open toe**

Putting on your compression stockings with an open toe is much easier with the included Ofa donning aid (100% polyester). Slip your foot into the dressing aid so that it reaches just above the heel.


Now pull the stockings on until the heel is in the right place. This is very easy with the Ofa dressing aid.

To take out the dressing aid, pull back the forefoot section of the stockings.

Pull the dressing aid out over the tip of your toes.

Then carefully roll back the forefoot section of the stocking again. Then pull the rest of the stocking up evenly without overstretching it. Finally, adjust the stockings on the leg until there are no wrinkles.

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**Ofa Bamberg GmbH**