PUSH MED ANKLE BRACE

INTRODUCTION

Push med products provide effective solutions for the treatment and prevention of injuries of the loconotor apparatus. Approved by physicians and scientists.

Your new ankle brace is a Push med product, characterised by a high level of effectiveness, comfort and quality. The Push med ankle brace increases ankle stability and helps control pain. This increases physical strength. Mechanical support, adjustable pressure and good position retention constitute the base of its performance.

Push med products are innovative and are manufactured for you with great care. This ensures maximum security and more freedom of movement. To promote satisfied use.

Read the product information carefully before you use the Push med ankle brace. If you have any questions, please consult your reseller or send an e-mail to info@push.eu.

PRODUCT CHARACTERISTICS

Your Push med ankle brace offers a high level of support and protection during the treatment or prevention of ankle injuries. It can be used in a wide range of cases and conditions. The combination of a zip and hook and loop fastening ensures that the sock can be fitted easily. Non-elastic straps can be used to set the pronation and supination limits, without any impediments for walking. The elastic straps with heel-lock, according to the functional tape concept, offer individually adjustable compression around the ankle. The silicone application on the inside of the bracket stimulates the ankle's proprioception. This is done to anticipate an imminent trauma.

The ankle brace has a perfect fit and can be worn in practically any shoe. Comfort zones of padded material have been created around vulnerable areas such as the instep and the heel.

The brace is made of Sympress™, a high-quality comfortable material. The use of microfibres makes the inside feel soft and keeps the skin dry. The Push med ankle brace is easy to apply and can be washed without any problem.

INDICATIONS

- Moderate to serious instability
- Acute lateral ankle ligament trauma
- Stable ankle fracture
- Follow-up treatment of ankle fractures
- Posturgical treatment of ligament injuries or fractures
- Capsular irritation
- Tendinitis
- Arthrosis with instability
- Secondary prevention of ankle ligament trauma
- Impingement (capsular clamping; anterior and posterior)
- Football ankle

There are no known disorders where the brace cannot or should not be worn. In the case of poor blood circulation in the limbs, a brace may cause hyperaemia. Avoid applying the brace too tightly in the beginning.

Pain is always a sign of warning. If the pain persists or the complaint deteriorates, it is advisable to consult your doctor or therapist.

DETERMINING SIZE / FITTING

The ankle brace is available in 5 sizes. Position the foot at an angle of 90° in relation to the lower leg. Then measure the circumference around the heel and the instep of the foot, as shown in the illustration on the packaging. The table on the packaging shows the correct brace size for your foot.

- A. Slide the base around your foot, making sure it fits snugly. Keep the foot at an angle of 90° in relation to the lower leg. Then close the zip and the elastic fastening at the front.
- B. Close the non-elastic straps crosswise across the instep. The correct fitting sequence is determined by the increasing molecule symbols on the hook and loop fasteners.
- C. The elastic straps are wound around the foot in the shape of the figure 8. Start by tensioning the heel-lock. Position the strap diagonally across the instep in the direction of the inner foot. Then pull the strap underneath the foot towards the outside and diagonally across the instep to the inside of the lower leg.
- D. Close the elastic strap with a single wrap around the lower leg.
- E. Tension the heel-lock of the second strap. Position the strap diagonally across the instep in the direction of the outer foot.
- F. Then pull the strap underneath the foot towards the inside and diagonally across the instep to the outside of the lower leg.
- Close the elastic strap with a single wrap around the lower leg.

WASHING INSTRUCTIONS

Before washing, close the hook and loop fasteners to protect the brace and other laundry. Close the brace can be washed at temperatures up to 30 °C, on a delicate fabrics cycle or by hand. Do not use any bleaching agents.

A gentle spin-dry cycle can be used; after this, hang the brace up to dry in the open air (without heating source or in a tumble dryer).

To extend the brace's service life, it should be dried after use, if necessary, before being stored.

IMPORTANT


- Proper use of the Push brace may require diagnosis by a doctor or therapist.
- Push braces were designed to be worn on an undamaged skin. Damaged skin must first be covered with suitable material.
- Optimal functioning can only be achieved by selecting the right size. Fit the brace before purchasing, to check its size. Your reseller's professional staff will explain the best way to fit your Push brace. We recommend that you keep this product information in a safe place, to be able to consult it for fitting instructions whenever needed.
- Consult a brace specialist if the product requires individual adaptation.
- Check your Push brace each time before use for signs of wear or ageing of components or seams. Optimal performance of your Push brace can only be guaranteed if the product is fully intact.
- Within the framework of EU regulation 2017/745 on medical devices, the user must report any serious incident involving the product to the manufacturer and to the competent authority of the Member State in which the user is established.

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