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Instructions for use

ofa[®]

**Produktinformationen und
Gebrauchsanweisung - online**
Product information and
instructions for use - online

dynamics[®]

Hallux Valgus Nachtschiene

REF 44 6005

- EN** Hallux Valgus Night Splint
NL Nachtsplank voor Hallux Valgus
IT Correttore Notturmo per Hallux Valgus
PT Tala Noturna para Hallux Valgus
SV Nattskena för Hallux Valgus
EL Νάρθηκας νυκτός βλαισού μεγάλου δακτύλου κάτω ποδός

1

2

deutsch	4
english	6
nederlands	9
italiano	11
português	14
svenska	16
ελληνικά	19

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english

Fitting instructions

The pictures that go with the fitting instructions can be found on the inside cover (2) at the start of the brochure.

1 Place your big toe in the front part of the splint so that the wide leather strap lies against the outside of the ball of your foot.

2 To secure the splint, guide the narrow Velcro strap around the outside of your foot and through the notch in the splint on the top of your foot. Tighten it until the required correction position is reached, and fasten the strap. Ensure that the product is secure but not so tight that it causes pinching.

Intended Use

The product is intended for the con-

servative correction of the metatarsophalangeal joint of the big toe in cases of hallux valgus, bunions and other deformities of the big toe as well as to provide stabilisation after hallux valgus surgery.

Important information

- The initial fitting and adjustment of the product as well as instructions on its proper and safe use should be provided by trained medical/orthopaedic specialists.
- The product should only be used for the indications in the adjacent column.
- Ensure that the product fits correctly and check it before each use for functionality, wear and damage (e.g. breaks, tears or deformations). A damaged product may no longer be worn.
- Wear the product directly on the skin. It should only be worn over intact skin. If it is to be worn in combination with other products, please discuss this with your doctor first.
- If you experience an unpleasant sensation, severe pain or other symptoms while wearing the product, remove it and contact your doctor or specialist supplier.

- We use a comprehensive quality assurance system to check our products. If you have any complaints, however, please contact your specialist supplier.

Instructions for care

- The product should be washed before its first use and then on a regular basis.
- The splint can be cleaned using a damp cloth.
- To maintain the quality of your product, please do not use any chemical cleaning products, bleaches, petrol or fabric softener on it. Do not apply any fatty or acidic creams, ointments or lotions to the skin in the area where the product is used. These substances may damage the material.
- Store the product in a dry place and protect it from sunlight and heat, preferably in the original packaging.

Care

-  Hand wash
-  Do not bleach
-  Do not tumble dry
-  Do not iron
-  Do not dry clean

Do not use fabric softener!

Reporting obligation

In accordance with legal provisions within the EU, patients and users are obliged to report any serious incident that occurs when using a medical device to both the manufacturer and the competent national authorities (in Germany BfArM, Bundesinstitut für Arzneimittel und Medizinprodukte, the Federal Institute for Drugs and Medical Devices) immediately.

Indications

- Hallux Valgus (bunion)
 - Formation of bony lumps
 - Deformed big toe
 - For stabilisation after bunion surgery
-

Material composition

38% Acrylonitrile butadiene styrene
22% Polyurethane (foam)
20% Polyamide microfibre
18% Polyamide
2% Polyurethane

Side effects

There are currently no known side effects for any part of the body if the product is used correctly. However, if the product is too tight it can cause local pressure points or constrict

blood vessels and nerves.

Contraindications

Consult your doctor before using this product in the following cases:

- Skin diseases or injuries in the area where the product is used, particularly if you experience any signs of inflammation such as redness, an increase in temperature or swelling
 - Feeling of numbness and circulatory problems in the area where the product is used
 - Problems with lymphatic drainage and unexplained swelling
-

Disposal

Please dispose of the product according to local regulations after the end of use.

Guarantee

If used as intended and maintained properly, the product can be used for six months. The product should be handled carefully in line with these instructions for use. Improper use or improper changes to the product may negatively impact the performance, safety and functionality of the product and lead to the exclusion of claims. The product is designed for the treat-

ment of a single patient. It is not designed to be reused by a second patient or any further patients.

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Ofa Bamberg GmbH