---

memory<sup>®</sup>

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EN

**Memory – the elegant medical compression**

Memory compression stockings shine with a transparent fine stocking look and an attractive colour palette. Ten standard colours and four seasonal colours make the collection a fashionable experience for elegant and sophisticated women. Memory is thus the perfect accessory for all occasions: whether with modern everyday clothes, functional business attire, or elegant evening wear.

**Memory Aloe Vera – the elegant medical compression, including skin care**

Memory Aloe Vera combines a transparent fine stocking look and elegant appearance with integrated skin care. Aloe vera extracts and vitamin E naturally moisturise the skin every time it is worn. Functional yarns and a sophisticated knitting technique provide a pleasantly soft feel. Memory Aloe Vera: beautiful, healthy and well-groomed legs all day long.

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---

EN

**Lastofa – the versatile medical compression**

The comfortable all-rounder wear in round-knitted: opaque and modern, with an exact fit and soft feel! Lastofa offers reliable support for men and women. In many variations and with thoughtful details. Also suitable for soft connective tissue and advanced venous insufficiency. In compression class 1 and 2 as standard and made-to-measure products.

**Lastofa with cotton – the versatile medical compression with cotton**

The robust and skin-friendly compression stocking with cotton directly on the skin. Especially suitable for patients who suffer from allergies and patients who prefer a robust quality. In compression class 2 and 3 as standard and made-to-measure products.

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# english

## Dear User,

we are delighted that you have chosen Ofa Bamberg compression stockings and thank you for putting your trust in us. Please take a moment to read through this product information carefully. For optimal compression therapy effectiveness!

## Intended use

Compression stockings are used to treat phlebological and lymphological arm and leg conditions, as well as lipoedema. Different products are used depending on the indication.

## Effective principle of medical compression stockings

Products tailored to the respective indication and controlled pressure accelerate blood reflux, improve microcirculation in the tissue and increase lymph flow.

## ⚠ Important information

Only wear compression stockings after consultation with your doctor. Wearing your stockings daily and holding regular consultations with your doctor will ensure successful treatment! The stockings must be obtained from a medically trained specialist retailer. Your specialist retailer will also advise you on how to put the compression stockings on correctly. If you experience leg pain or skin irritation while wearing these stockings, please consult your doctor immediately. You can prevent skin reactions in the area around the grip tops by slightly adjusting the grip tops several times during the day. Please contact your doctor about wearing the stockings in combination with using other medical devices. We test our products within the scope of a comprehensive quality assurance. However, if you still have cause for complaint, please contact your specialist retailer. Do not attempt to perform repairs yourself. Return the stockings to the specialist retailer who sold them to you. Any pilling of the material that may arise due to friction with other garments will not reduce the stocking's quality.

## Care instructions

The product is suitable for immediate use. We recommend washing the product before wearing it for the first time. Then wash daily by hand or on a delicate fabric program. Wash the stockings separately with mild detergent without fabric softener (e.g. Ofa Clean special detergent) at a maximum temperature of 40 °C and rinse them well. The stockings can be spun gently, but do not wring them out. Dry your compression stockings in air or in a dryer on a gentle cycle. Do not place them on a heater or in direct sunlight. Never iron the stockings. To maintain the effectiveness of your stockings, please do not use chemical cleaning agents, petrol, or fabric softeners. These substances damage the knitted fabric. Compression stockings are largely resistant to fats and oils, but should not come into contact with ointments. Stockings with silicone-coated grip tops should generally not come into contact with solvents, ointments, re-greasing washing lotions, and body creams.

## Indications (see also AWMF Medical Compression Therapy Guideline 037/005, Version 12/18)

- Improvement of venous symptoms and quality of life with chronic venous diseases
- Prevention and treatment of venous oedema
- Prevention and treatment of venous skin alterations
- Eczema and pigmentation
- Lipodermatosclerosis and atrophic blanche
- Therapy and pain reduction for venous leg ulcers, including relapse prevention
- Treatment of arterial and venous leg ulcers (see contraindications!)
- Varicose veins
- Initial phase after treatment for varicose veins
- Functional venous insufficiency
- Venous malformations
- Venous thrombosis
- Condition after thrombosis
- Post-thrombotic syndrome
- Thrombosis prophylaxis in mobile patients
- Lymphoedema
- Oedema in pregnancy
- Post-traumatic/post-operative oedema
- Cyclic idiopathic oedema
- Stage II lipedema and above
- Venous stasis resulting from immobility (arthrogenic congestive syndrome, paresis, and partial paresis of the extremity)
- Occupation-related oedema (occupations involving long periods of standing and sitting)

- Medication-related oedema, if change of medication not possible
- Obesity with functional venous insufficiency
- Inflammatory dermatoses of the legs
- Nausea, dizziness, and congestion complaints during pregnancy
- Condition after burns
- Scar treatment

## Contraindications

### Absolute contraindications:

- Advanced peripheral arterial occlusive disease (if any of these parameters apply: ABPI < 0.5, ankle artery pressure < 60 mmHg, toe pressure < 30 mmHg, or TcPO<sub>2</sub> < 20 mmHg arch)
- Congestive heart failure
- Septic phlebitis
- Phlegmasia cerulea dolens

### Risks to be considered

- (therapy decision based on benefit/risk assessment)
- Weeping dermatoses
- Allergies to compression stocking materials
- Severe sensitivity disorders of the extremity
- Advanced peripheral neuropathy (e.g. diabetes mellitus)
- Chronic polyarthritis

## No liability is accepted if these contraindications are not observed.

## ⚠ Risks and side effects

Medical compression stockings can cause skin necrosis and pressure damage to peripheral nerves. Patients with sensitive skin may experience itching and skin irritation when using compression stockings. Adequate skin care is therefore recommended.

If you suffer the following symptoms, please take off the stockings immediately and discuss further treatment with your doctor: blue or white coloration of the toes, sensations of discomfort and numbness, increasing pain, shortness of breath and sweating, acute movement restrictions.

## Warranty

The shelf life of medical compression stockings is limited to 3 years when properly stored. Added to this is the wearing time, which is a maximum of 6 months with regular use and proper care. After this period, the products are excluded from the warranty. Defects attributable to normal wear and tear or improper use are not covered by the guarantee. Improper use also includes any subsequent modification to the product that has not been carried out by the manufacturer.

The product is intended for use by a single patient. If it is passed on to other persons for reuse, the manufacturer's product liability is no longer valid.

## Disposal

After the end of the usage period, please dispose of the product in accordance with local regulations.

## Mandatory reporting

Due to legal requirements, patients and users are obliged to report immediately any serious incidents which have occurred in connection with the product and which have resulted in, or could have led to a significant deterioration in the state of health or to death, to the manufacturer and the competent national authority (in Germany BfArM).

## english

Otherwise, the adhesive strength of the grip tops could be permanently impaired (leg hair and dry skin can also limit the function of the grip tops). Store the stockings in a dry place away from sun and heat, in the original packaging if possible.

### **Wash inside out and do not use fabric softener!**

- Wash on a delicate or gentle cycle
- Do not bleach
- Dry gently and carefully
- Do not iron
- No chemical dry cleaning

### **Material composition**

Polyamide (Polyamid), elastane (Elastan) and, for some products, cotton (Baumwolle) is used. Please refer to the textile label sewn into your stocking for the precise details.

### **Tips for donning the stockings**

Put on the stockings first thing in the morning after getting up, as long as your legs are free from swelling. Please do not pull the compression stocking up by the edge or the grip top, instead, slide it upwards with your hands. We recommend wearing household rubber gloves to make it easier to put on the stockings and to protect the knitted fabric. Jewellery and long fingernails can damage the stockings. Please also pay attention to regular foot care and proper footwear. Ask your medical retailer about Ofa Grip special gloves and Ofa Fit Magnide, the practical donning and doffing aid for compression stockings from Ofa Bamberg, available with open and closed toe. Ofa Fit Magnide is made of very smooth canvas material and enables the stockings to slide easily over the foot.

### **How to put them on, closed toe**

To put the stockings on, turn it inside out as far as the foot. Slip the foot section over your toes and then over the instep.

Now carefully pull the leg section over your heel. Please avoid pulling up the stockings at the end or grip top.

Pull the rest of the stockings up evenly without overstretching it. Finally, adjust the stockings on the leg until there are no wrinkles.

### **How to put them on, open toe**

Putting on your open toe compression stockings is much easier with the supplied Ofa dressing aid. Slip your foot into the dressing aid so that it reaches just above the heel.

Now pull the stockings on until the heel is in the right place. This is very easy with the Ofa dressing aid.

To take out the dressing aid, pull back the forefoot section of the stockings.

Pull the dressing aid out over the tip of your toes.

Then carefully roll back the forefoot section of the stockings again. Then pull the rest of the stocking up evenly without overstretching it. Finally, adjust the stockings on the leg until there are no wrinkles.

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Ofa Bamberg GmbH