## LIGADURA DE PULSO PUSH SPORTS

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## ENGLISH (EN)

The Push Sports Wrist Brace stabilises the wrist joint and keeps the palm of the hand free to hold or make contact with any piece of sports equipment.

The brace restricts movement, without immobilising the joint, allowing you to continue playing sports. The elastic strap enables you to determine the degree of compression that's right for you. The wrist brace is easy to apply with one hand, and the quick-drying material maximizes comfort. If you have any questions, please consult your retailer or send an e-mail to [info@push.eu](mailto:info@push.eu).

### WHEN TO USE THE PUSH SPORTS WRIST BRACE?

Wrist pain may be caused by overstraining the tendons in the wrist and hand. Particularly effective in the event of recurrent wrist pain, the Push Sports Wrist Brace decreases the strain on the active but vulnerable wrist.

#### INDICATIONS

- To prevent recurring strain symptoms in the wrist joint
- Strain symptoms in the wrist joint
- Pain in the wrist joint
- Laxity or instability in the wrist joint

#### CONTRAINDICATIONS

- There are no known disorders where the Push Sports Wrist Brace should not be worn.

#### CAUTIONS

- Ensure that the brace is not fitted too tightly. Use the brace only for the appropriate indications. Proper use of your Push Sports brace may require diagnosis by a doctor or therapist. Pain is always a warning sign. If pain or complaints persist without a clear cause, it is recommended that you visit your GP, physiotherapist or sports physician.
- Push Sports braces were designed to be worn on an undamaged skin. Damaged skin must first be covered with suitable material.
- Optimal functioning can only be achieved by selecting the right size. Therefore, it is recommended that you fit the brace before purchasing to check its size. When buying online, always ensure that there is a good return service. Your reseller's professional staff will explain the best way to fit your Push Sports brace. We recommend that you keep this product information in a safe place, to be able to consult it for fitting instructions whenever needed.

- Check your Push Sports brace each time before use for signs of wear or ageing of components or seams. Optimal performance of your Push Sports brace can only be guaranteed if the product is fully intact.
- Within the framework of EU regulation 2017/745 on medical devices, the user must report any serious incident involving the product to the manufacturer and to the competent authority of the Member State in which the user is established.

#### SPECIAL FEATURES

The Push Sports Wrist Brace provides relief by decreasing the strain on the active but vulnerable wrist. The wrist brace is fitted with stiffeners that restrict the bending of the wrist joint. The brace is closed with an elastic strap that keeps the stiffeners in place, provides comfortable compression around the wrist joint, and reinforces the supportive action of the stiffeners.

#### TAKING MEASUREMENTS

The Push Sports Wrist Brace is available in three sizes, and in right and left versions. Select the correct brace size by measuring the circumference of your wrist.

#### FITTING

- ① Insert the thumb in the opening provided. The D-ring goes at the back of the hand.
- ② Tighten the elastic strap around the wrist joint and fasten the hook and loop fastener. The buckle makes it easy to tighten the elastic strap.

[pushsports.eu/wrist-brace/applying](https://pushsports.eu/wrist-brace/applying)

For details of the distributor in your country, please refer to [www.push.eu/distributors](http://www.push.eu/distributors)

#### WASHING INSTRUCTIONS/MAINTENANCE

To promote hygiene and to maximise the lifespan of the brace, it is important to wash the brace regularly, observing the washing instructions. The brace can be washed at temperatures of up to thirty degrees Celsius, on a delicate fabric cycle or by hand. Before washing, fasten the hook and loop material to protect the brace and other laundry. Do not use any bleaching agents. Let the brace dry in the open air (not near a source of heat or in a tumble drier). Never store the brace when it is wet. Always let it dry first.

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